

GREAT PLAINS OAT



Please read the instructions thoroughly provided with this kit. Check your kit and call our office or the lab if anything is missing or damaged.

Instructions

- Freeze the cold pack the night before collection.
- Collect sample from first urine of the morning prior to eating or drinking. Begin capture midstream to ensure good quality sample. Once collected, label the sample container with name, date of birth, date and time of collection, and place in bio-hazard bag to ship.
- Please be sure the requisition form and the accompanying paperwork are filled out before sending and include a copy of your Medicare card if applicable.

Ensure that the requisition and all accompanying paperwork are filled out and placed in the package before sealing.

