

SANESCO NEUROTRANSMITTERS



Please read the instructions thoroughly provided with this kit. Check your kit and call our office or the lab if anything is missing or damaged.

- Urine Collection - DO NOT collect the first urine of the morning.
- Saliva Collection – Collect first saliva sample one hour after rising. Samples are typically collected before breakfast, before lunch, before dinner and before bed.
- Please ensure that the completed requisition form and symptom questionnaire are included before sealing to ship and include a copy of your Medicare card if applicable.

Ensure that the requisition and all accompanying paperwork are filled out and placed in the package before sealing.

